

**Title Page.** *Required.*

**Abstract.** *Required.*

**Introduction.** *Required.*

Please address the following with regard to the lab:

1. Describe some applications of the titration that you will be performing.
2. Describe the electrochemistry relevant to the calculations. Be sure to include the following concepts/expressions:
  - a. The shorthand notation for the electrochemical cell at the equivalence point
  - b. The balanced electrochemical reaction corresponding to the cell
  - c. The expression for the overall standard electrochemical potential
  - d. The Nernst equation describing the electrochemical potential
  - e. The chemical equilibrium expression between  $\text{AgCl}$ ,  $\text{Ag}^+$ , and  $\text{Cl}^-$  before the equivalence point when there is an excess of  $\text{Ag}^+$
  - f. The expression for  $K_{\text{sp}}$
  - g. The expression for the cell potential in terms of  $K_{\text{sp}}$ , 0.100 M  $\text{CuSO}_4$ , and the overall standard electrochemical potential
3. Include a labeled picture of the entire set up and show the electron flow. Directly address points 2a-d.

**Experimental.** *Required.*

**Results.** *Required.*

1. Include titration curves and a first derivative plot of all titrations performed.
2. Include tables containing key pieces of the raw data (volume needed to reach the equivalence point,  $E_{\text{cell}}$  at the equivalence point, etc.)
3. Calculate the concentration of the silver nitrate solution using the data from the known concentration of NaCl.
4. Calculate the  $K_{\text{sp}}$  of AgCl two ways:
  - a. Using the potential at the equivalence point of the titration of the **known concentration** of NaCl solution.
    - i. Use five other students' data so that a 95% CI can be calculated.
    - ii. Compare to known value of  $K_{\text{sp}}$  ( $1.8 \times 10^{-10}$ ) by calculating a relative error.
  - b. Using the **unknown concentration** of NaCl solution.
    - i. Use five other students' data so that a 95% CI can be calculated.
5. Calculate the concentration of chloride in the unknown two ways:
  - a. Using the titration curve and the volume used to reach the equivalence point (similar calculation to an acid/base titration).
    - i. Use five other students' data so that a 95% CI can be calculated.
  - b. Using the initial  $E_{\text{cell}}$  and the known value of  $K_{\text{sp}}$ .
    - i. Use five other students' data so that a 95% CI can be calculated.

**Discussion.** *Required.*

**Conclusion.** *Required.*

## **Lab Report Requirements, Week 6**

*Revised: 5/21/15*

**References.** *Always required.* Lab manual website must be cited.  
**Statistical Analysis.** *Required.*